NISSEI DSK-1031

INSTRUCTOIN MANUAL DIGITAL BLOOD PRESSURE MONITOR DSK-1031

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
INDICATIONS FOR USE

This Instruction Manual is designed to assist the user with safe and effective operation of the digital blood pressure monitor DSK-1031 (hereinafter – the “Device”). The device is designed to measure systolic and diastolic blood pressure and heart rate reading in patients aged 12 years and older. This device should not be used for newborn children. The incorrect measurement result is also possible when the device is operated in a baby hand. Consult your doctor for blood pressure measurement in children. The device is recommended for use in patients with unstable (nonpermanent) blood pressure or hypertension at home as a supplement to medical surveillance. The cuff is suitable for an upper arm with a circumference of about 22 - 42 cm. Blood pressure is measured in the range from 40 to 250 mm Hg and heart rate, in the range from 40 to 160 heartbeats per minute.

OPERATION PRINCIPLES

The device uses the oscillometric method of measurement. The cuff connected to the electronic module is wrapped around the upper arm.

When you press the START/STOP button, the device starts the automatic inflation, which is followed by the blood pressure measurement. The sensor element of the device detects slight pressure oscillations in the cuff produced by the expansion and contraction of the brachial artery in response to each heartbeat. Inflation is stopped when the cuff is pumped enough to determine diastolic and systolic pressure, whereupon the air is discharged from the cuff. Rhythm and amplitude of the pressure waves are measured and displayed on the LC display as a numerical value in millimeters of mercury. The device has an arrhythmia indicator, as well as two memories with 60 cells in each calculating the average value.

NISSEI New Technologies

Fuzzy Inflation Algorithm – is an algorithm for automatic selection of the cuff inflation pressure. Using this algorithm, the device by itself determines the pressure level to which it is necessary to inflate the cuff based on the patient’s systolic pressure. Owing to the Fuzzy Inflation algorithm the device becomes easier to use, while the measurement gets more comfortable and more accurate.

Arrhythmia indicator is a special icon on the display that informs on the irregular heartbeat, while the measurement result is correct.

Touch control – is control of the device by slight touching it.

Detection of interference noises – indicator informs on the occurrence of noises that can affect the measurement result.

Control of correct fixing the cuff – shows if the cuff is fastened too tight or too loose.

Indication of reliability – this symbol is displayed when all stages of proper measurement procedure are observed.

Pulse pressure – along with the measurement result the device displays the pulse pressure value. Pulse pressure is a difference between systolic and diastolic pressure.

WARNING! Using a cuff that is different from that supplied with this device is not allowed.

COMPLETE SET

The complete set includes

- Electronic unit - 1 pc.
- Cuff (including air hose and air hose plug) - 1 pc.
- Batteries - 4 pcs.
- Power supply source - 1 pc.
- Bag - 1 pc.
- Instruction Manual - 1 pc.
- Warranty - 1 pc.
- Packaging - 1 pc.

RECOMMENDATIONS ON CORRECT MEASUREMENTS

1. If treated with hemodialysis or anticoagulants, antiplatelets or steroids, refer to your doctor about the blood pressure measurement.
2. Malfunctions are possible when the device is used near working mobile phones, microwave ovens and other equipment generating electromagnetic radiation.
3. For correct measurement it is necessary to know that the BLOOD PRESSURE IS SUBJECT TO SHARP FLUCTUATIONS EVEN IN SHORT TIME INTERVALS. The blood pressure level depends on many factors. It is commonly lower in summer and higher in winter. Blood pressure varies along with atmospheric pressure and depends on the physical exertion, emotional excitability, stress and diet. Medical drugs, alcohol and smoking exert great influence as well. Occasionally, measurements in the clinic cause

an increase in pressure values. Therefore, blood pressure measured at home is often different from that measured in the clinic. Since blood pressure increases at low temperatures, measurements should be made at room temperature (about 20°C). If the device was stored at low temperature, prior to using it should be kept for at least 1 hour at room temperature, otherwise the measurement result can be erroneous. During the day, the difference in the readings in healthy people may attain 30-50 mm Hg for systolic (upper) pressure and up to 10 mm Hg for diastolic (lower) pressure. Dependence of blood pressure on various factors is individual for each person. Therefore it is recommended to keep a special recording of blood pressure readings. **ONLY A DOCTOR MAY ANALYZE TRENDS IN CHANGING YOUR BLOOD PRESSURE BASED ON CORRESPONDING RECORDINGS.**

4. In case of cardiovascular diseases and a number of other diseases that require the blood pressure monitoring, measurements should be carried out in the hours specified by a doctor. **REMEMBER THAT THE DIAGNOSTICS AND ANY TREATMENT OF ARTERIAL HYPERTENSION SHOULD BE CARRIED OUT ONLY BY A DOCTOR BASED ON BLOOD PRESSURE READINGS OBTAINED BY A DOCTOR. MEDICAL DRUG ADMINISTRATION OR CHANGE OF DOSAGES SHOULD BE MADE ONLY BY PRESCRIPTION OF AN ATTENDING DOCTOR.**


Fig.1

5. In case of disorders such as deep vascular sclerosis, weak pulse wave and break in rhythm of heart contractions, the correct blood pressure measurement can be complicated. **IN THIS CASE, A DOCTOR SHALL PROVIDE RECOMMENDATIONS IN RELATION TO USE OF THIS DEVICE.**
6. **KEEP QUIET DURING THE MEASUREMENT TO OBTAIN THE CORRECT BLOOD PRESSURE READING WHEN USING THE ELECTRONIC DEVICE.** The blood pressure meas-

urement should be carried out in a quiet comfortable atmosphere at room temperature. Exclude meal an hour before the measurement, and exclude smoking, soft drinks, and alcohol 1.5-2 hours before the measurement.

7. Accuracy of the blood pressure measurement depends on matching the device cuff and size of your hand. THE CUFF SHOULD NOT BE TOO SMALL OR TOO BIG.

8. Repeated measurements are carried out at 5-minute intervals to recover the blood circulation. However, persons suffering from severe atherosclerosis, due to a significant loss of elasticity of blood vessels, need longer intervals between measurements (10-15 minutes).

This also concerns patients suffering from long-term diabetes. For more accurate determination of blood pressure it is recommended to carry out a series of three consecutive measurements and to calculate the average value of measurement results.

POWER SUPPLY OF THE DEVICE

INSTALLATION OF BATTERIES


Fig.2


Fig.3

1. Open the battery compartment (Fig. 2).

2. Install four "AA" batteries in the compartment.

Make sure that polarity corresponds to signs (+) and (-) shown inside the compartment (Fig. 3). Batteries are readily installed by pressing the end "-" on the spring.

Allowed to rechargeable batteries. To charge the batteries, use the special charger (not included).

3. Close the battery compartment. Do not use excessive force when removing the cover.

Do not use excessive force when removing the cover.

Battery Replacement Indicator

Replace all the batteries when the battery replacement indicator is flashing on the display during the measurement. If upon the device turning on the indicator is steadily flashing, the measurement will not be possible until all the batteries are replaced. The battery replacement indicator does not show a discharge degree.

Use alkaline batteries to increase the device operation duration. Ordinary zinc-carbon batteries require more frequent replacement. The enclosed batteries are meant for testing the sold device, and their operation period can be less than that of batteries acquired in the trade network.

Since neither the device nor the batteries are the waste that can be utilized at home, follow your national/local regulations for waste recycling and take them to corresponding collection facilities.

USE OF THE DEVICE WITH THE POWER SOURCE

Socket for the power source is arranged on the left side of the device. To use the device with the power source, connect the power source connector to the device, install the power source plug into the socket outlet, and press the «START/STOP» button.

When finished, turn off the device by pressing the «START/STOP» button, unplug the power source from the socket outlet and disconnect the power source connector from the device.

NOTE!

If there is no battery in the device, turning off the power source will result in zeroing of measurement results stored in the device memory and set date and time. If you do not want to make the data erased, do not remove the batteries from the device when using the power source.

SETTING DATE AND TIME

Date and time can be set after installing batteries. Setting the date and time guarantees the preservation of measurement results with indicated correct date and time. The device can be used without setting the date and time.

Press and hold the SET button until the display flashes the value of the year.

Date and time are set in the following order: year, month, day, hour and minute.

1. Setting the Year

Use the **[1]** button to increase and **[2]** button to decrease the year value. Press the SET button to confirm and to pass to the next step.

2. Setting the Month

Use the **[1]** button to increase and **[2]** button to decrease the month value. Press the SET button to confirm and to pass to the next step.

3. Setting the Day

Use the **[1]** button to increase and **[2]** button to decrease the day value. Press the SET button to confirm and to pass to the next step.

4. Setting the Time

Watch uses a 12-hour time format of day. Use the **[1]** button to increase and **[2]** button to decrease the hour or minute value. Press the SET button to confirm the settings. To stop the setting, press the «START/STOP» button.

IMPORTANT! If the date and time are set, then when turned off the device display will show the current time.

CORRECT POSITION DURING MEASUREMENT

Sit down at the table so that during the blood pressure measurement your hand is on its surface. Make sure that the place where the cuff is put on the upper arm is about the same level as the heart and the hand freely lies on the table and does not move (Fig. 4).


Fig.4

Fig.5


Fig.6

You can also measure your blood pressure when lying on your back. Look up, stay calm and do not move during the measurement. Make sure that the place where the cuff is put on the upper arm is about the same level as the heart (Fig. 6).

Measured values may vary slightly, depending on the position during the measurement. If the cuff is above/below the level of the heart, resulting reading may be incorrect (lower/higher)

CUFF PREPARATION


Fig.7


Fig.8


Fig.9


Fig.10


Fig.11

1. Apply the cuff to your left upper arm so that the air tube is directed to your palm (Fig. 7). If the measurement on your left arm is difficult, you may use your right arm. In this case remember that the readings may differ by 5-10 mmHg and even more.

2. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. Air tube should be directed towards the palm (Fig. 8).

3. Fix the cuff so that it fits tightly to the arm, but see that it is not overtight (Fig. 9). Too tight or too free placement of the cuff may give inaccurate readings.

4. If the hand is cone-shaped, it is recommended to put the cuff spirally, as shown in the figure (Fig. 10).

5. If the rolled-up sleeve squeezes the arm interfering with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure (Fig. 11).

MEASUREMENT PROCEDURE

IMPORTANT! The device has touch-sensitive buttons that are pressed by slight touching. Moisture, dirt and extraneous objects between the finger and the device panel can affect the ability of buttons to respond to touching.

1. Insert the Air Tube Plug into the Cuff Jack.

Before measurement, take a few breaths and relax. During the measurement, do not talk and do not move.

2. Press the START/STOP button. The deflating symbol «V» will flash on the display, and the device will release the remaining air from the cuff (Fig. 12).

3. Beep will be heard, and air will be pumped into the cuff. The symbol «^» will flash, and the displayed value will increase (Fig. 13). Inflation will stop at optimum level owing to the Fuzzy Inflation algorithm.

4. The symbol «^» will disappear, and the measurement will start. The cuff pressure will decrease slowly.


Fig.12


Fig.13

Noise Interference Detection

This device detects the presence of extraneous sounds or noises that can affect the measurement results. When such interference is detected, the symbol «⚡», appears on the display. In this case, it is recommended to repeat the measurement to eliminate the interference.

Control of Correct Fixing the Cuff

The symbol «OK» indicates that the cuff is applied correctly, the symbol «⊖» indicates that the cuff is fastened too tight or too loose. In this case, it is recommended to repeat the measurement after proper fixing the cuff.


Fig.14

Press the START/STOP button to stop forcibly the measurement: the device will stop inflation and quickly release the air.

Automatic Additional Inflation.

If on starting the measurement, the operating pressure appears to be insufficient due to hand movements or due to excessive tension in a hand, the cuff will be filled with air again to the level of about 30 mm Hg higher than the original level. Automatic additional inflation will be repeated until the measurement is completed successfully. This is not a malfunction.

5. The symbol «♥» will start flashing synchronously with a sound signal, as soon as the sensor catches the pulse pressure (Fig. 14).

6. When the measurement is completed, the display will show the values of blood pressure, pulse pressure, WHO scale and heart rate (Fig. 15). The device will automatically release the air from the cuff.

7. Press the **1** button or **2** button and the result will remain in the selected memory unit.

8. Press the START/STOP button to turn off the device.

If you forget to turn off the device, it will do so automatically after 3 minutes.

Pulse, stroke/min.

Systolic pressure, mmHg

Diastolic pressure, mmHg

The Cuff Symbol

The Reliability Symbol

Puls pressure, mmHg


Fig.15

Indication of Reliability

Important aspects in measuring the blood pressure imply the correct application of the cuff and proper posture during the measurement. These aspects may influence the measurement result. The device displays the symbol «☑», when the measurement procedure is carried out properly. If this symbol is not displayed, it is recommended to repeat the measurement observing the guidelines for proper measurement described in this manual.

Do not perform several measurements in a row.

This will cause numbing the hand and can affect the measurement result. Give your hand a break for at least 5 minutes.

INDICATION OF ARRHYTHMIA

The flashing symbol «⊕», appeared on the display reports on an irregular heartbeat rhythm (Fig. 16). If this symbol appears from time to time, refer to your doctor. The arrhythmia indicator occurrence can be also caused by body movements in the course of the measurement.


Fig.16

INDICATION OF READINGS BY THE WHO SCALE

In addition to pressure numerical values, the result is also displayed as a graphic scale. This scale makes it possible to classify the obtained blood pressure value in accordance with the World Health Organization recommendations. The scale appears along with the numerical value of blood pressure in the lower right corner of the display (Fig. 17).


Fig.17

Table of Results by the WHO Scale

Indication	WHO Classification	SYS	DIA
	Hypertension (severe)	≥180	≥110
	Hypertension (moderate)	160-179	100-109
	Hypertension (slight)	140-159	90-99
	Increased normal pressure	130-139	85-89
	Normal pressure	120-129	80-84
	Optimal	<120	<80

Indication of Pulse Pressure

This device calculates and displays the pulse pressure value (Fig. 18). Pulse pressure – a difference between the systolic and diastolic pressures – tends to increase with aging.

Systolic blood pressure continues to increase with aging, while diastolic blood pressure tends to decrease since the age of about 50 years. It is considered that a high pulse pressure is related to poor elasticity of arteries and is one of the risk factors for circulatory diseases. Pulse pressure is considered to be normal when it attains 35±10 mm Hg.


Fig.18

MEMORY FUNCTOIN

The measured values are automatically saved for later viewing in the selected memory unit. These memory units can be used to save separately the measurement results of two different persons or to save separately the results of morning and evening measurements.

Each unit can keep up to 60 measurement results and the average value thereof. When the number of kept values exceeds 60, the oldest readings are deleted to record new values.

Saved values are kept with the measurement date and time, if set. If you want to save time and date together with the measurement values, date and time should be set before the measurement.

When an error occurs (ERR), the measurement results are not saved.

Viewing the Saved Data

1. Turn off the device by pressing the START/STOP button. To view the results stored in the memory unit 1, press [1], and to view the results stored in the memory unit 2, press [2]. The selected memory unit will be underlined on the display.

The display shows the average value of stored results indicated with index «n» (Fig. 19). Average value is displayed when the selected memory unit contains two or more saved measurement results.

2. Each time you press the [1] button or [2] button, saved measurement results will be displayed in a sequential manner.


Fig.19

3. Memory unit number, measurement date and time will be displayed one by one in the upper right corner of the display.

4. The result saved in the memory unit 1 is the most recent among the saved data in the selected memory. The higher memory number, the older the result.

The memory data are displayed for about 30 seconds and after that, if no button is pressed, the device is automatically turned off.

Pressing the [2] button switches the display from records in the memory unit 1 to records in the memory unit 2, while pressing the [1] button results in returning to data recorded in the memory unit 1.

5. Press the START/STOP button to turn off the device.

Deletion of Saved Data



Readings can be deleted unit by unit: all together or separately. The memory may be cleared when the selected memory unit contains two or more results.

1. Select a value from the memory unit to be deleted or the average value (unit with index « $\bar{}$ ») to clear the whole memory unit.

2. Press and hold the [1] button or [2] until the sign «---» appears on the display.

INFORMATION ABOUT ERRORS

INDICATION	LIKELY CAUSE	METHODS OF CORRECTION
Blood pressure is too low or too high.	The cuff is not at the heart level. The cuff is put on incorrectly. During the measurement, a person was talking or moving.	Put on the cuff at the heart level. Check the cuff position on the hand. Be calm and quiet during the measurement.
Measurement results are different each time.	Effect of measurement conditions, physical or mental state.	Take measurements under the same conditions.
Measurement results are different in clinic and at home.	Effect of relaxed state at home and tension in clinic.	Show the pressure records made at home to your doctor for advice.
Inflation is repeated.	Pressure can not be measured due to movement or talking.	Re-inflation is not a malfunction of the device. Do not talk and do not move during the measurement.

	Maximum allowable pressure: pressure can not be measured due to movement or talking during the measurement, in spite of the fact that the cuff was pumped at maximum.	Do not talk and do not move during the measurement.
	Pressure can not be measured due to movement or talking.	Do not talk and do not move during the measurement.
	Cuff is not securely connected to the device. Cuff is not put on properly.	Check the connection. Make sure that the cuff is put on correctly.
No time indication on the display.	Time was not set. Note: there is no time indication when batteries and/or power source are not provided.	Set time and date. Install batteries and connect the power source.
Date and time are displayed as «-/-».	Time was not set or measurements were carried out before setting the time.	Set date and time. Date and time can not be saved without setting.
	Batteries are discharged.	Replace all batteries with new ones.
Nothing is shown on the display.	Batteries are discharged. Batteries are installed incorrectly. Connecting terminals are contaminated. Power source is not connected. Insufficient contact when touching the button.	Replace all batteries with new ones. Install batteries properly. Wipe connecting terminals with a dry cloth. Connect the power source. Ensure stronger pressing the button.

	You touched the START/STOP button when installing the batteries.	Turn off the device by pressing the START/STOP button and perform the measurement again.
---	--	--

If, despite the above-given recommendations, you fail to obtain the right measurements, stop the operation and contact the service center (addresses and telephone numbers of authorized organizations are provided in the warranty certificate). Do not attempt to adjust the device internal mechanism on your own.

WARRANTY

1. The manufacturer guarantees the warranty period of 5 years for the device from the sale date provided that the consumer observes operation, transportation and storage requirements. The warranty period for the cuff and the power source is 12 months from the sale date.
2. Warranty liabilities are documented with the warranty certificate upon selling the device to the buyer. The guarantee is valid provided that the device has not been opened or damaged by the buyer.
3. Addresses of organizations engaged in the warranty service are specified in the warranty certificate.

TECHNICAL SPECIFICATIONS

Operating Principle	Oscillometric method
Indicator	15 digits liquid crystal display
Pressure Indicating Range , mmHg	3-300
Measuring Range :	
cuff pressure, mmHg	40-250
pulse rate, bpm	40-160
Accuracy:	
cuff pressure, mmHg	±3
pulse rate, %	±5
Inflation	Automatic inflation (air pump, Fuzzy Inflation algorithm)
Deflation	Automatic Deflation (electric valve)
Electric supply voltage, V	6

Electric supply type	Four 1.5 volt LR6 (AA alkaline) batteries or AC adaptor
Power Consumption , W	4
Memory	2 units, 60 values per each + average value
ADP-W5 Power Source	
Output voltage, V	6
Maximum load current, A	0,5
Input voltage, V/Hz	100-240/50
Operating Condition	
temperature, °C	+10 to +40
humidity, % Rh	15-85
atmospheric pressure, hPa	700-1060
Transportation / Storage Condition	
temperature, °C	-20 to +60
humidity, % Rh	10-95
Cuff Size	adult (circumference of the upper arm 22-42 cm)
Overall dimensions:	
Size (electronic unit), mm	115 x 115 x 65
Weight (without packaging, bags, batteries and power source), g.	406
Manufacture year:	is specified on the device housing (in the battery compartment) in the device serial number after the letters «SN»
Protection level	IPX0 (IEC 60601-1)
Key to Symbols	<input checked="" type="checkbox"/> Type BF applied part <input checked="" type="checkbox"/> Refer to instruction manual/booklet <input checked="" type="checkbox"/> When utilizing the waste, refer to current rules applicable in your region

CARE, STORAGE, REPAIR AND DISPOSAL

1. This device should be protected from excessive moisture, direct sunlight, strokes, and vibration. **THE DEVICE IS NOT WATERPROOF.**
2. Do not keep or do not use the device in close proximity to heaters and open flame.
3. If the Device was stored at a temperature below the freezing point, keep it at least for 1 hour in some warm place before use.
4. If the device has not been used for a long time, remove the batteries. Leaking of batteries can cause damage to the device and terminate the warranty. **KEEP BATTERIES AWAY FROM CHILDREN!**
5. Do not contaminate the device and protect it from dust. The device can be cleaned with a dry, soft cloth.
6. Do not allow the contact between the device and its parts with water, solvents, alcohol, and gasoline.
7. Keep the cuff away from sharp objects, and do not try to pull out the cuff.
8. Do not expose the device to strong strokes and do not throw it.
9. The device does not contain any adjustment controls for settings. Unauthorized opening of the electronic device is forbidden. If needed, repair the device only in specialized organizations.
10. On the expiry of the specified operation term, refer to specialists (specialized repair organizations) on a periodic basis to check the technical condition of the device.
11. When utilizing the waste, refer to current rules applicable in your region. No special utilization conditions are specified by the manufacturer for this device.
12. The cuff is resistant to repeated sanitation. The cuff internal fabric surface (being in contact with a patient' hand) can be treated with a cotton swab moistened in a 3% solution of hydrogen peroxide. Partial discoloration of the cuff covering tissue is possible if used for a long time. Do not wash the cuff and do not treat it with a hot iron.
13. Do not leave unattended the device plugged into the network.

CERTIFICATION AND STATE REGISTRATION

The production of devices is certified pursuant to international standards such as ISO 9001, ISO 13485, ISO 14001:2004.

The device meets international standards IEC60601-1:2005+CORR.1:2006+CORR.2:2007, IEC60601-1-2:2001+A1:2004/IEC60601-1-2:2007 CISPR 11:2009+A1:2010.

✉ Complaints and requests should be addressed to:

Little Doctor Europe Sp. z o.o.

57G Zawila Street, 30-390, Krakow, Poland

Service phone: +48 12 2684748, 2684749

Produced by Nihon Seimitsu Sokki Co., Ltd.


Address: 2508-13 Nakago Shibukawa Gunma 377-0293 Japan

Manufacturer: Pt. Nss Indonesia

Address of Manufacturer: Blok A-2 No. 29 St4a Kawasan, Berikat Besland Pertiwi, Kota Bukit Indah Purwakarta 41181, Indonesia

Exporter: Little Doctor International (S) Pte. Ltd.

Address of Exporter: 35 Selegie Road #09-02 Parkline Shopping Centre, Singapore 188307

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Exportator: Little Doctor International (S) Pte. Ltd.

Adresa exportatorului: 35 Selegie Road #09-02 Parkline Shopping Centre, Singapore 188307

TECHNICAL DESCRIPTION

DSK-1031 complies with the EMC, electromagnetic compatibility, standard, IEC60601-1-2. Refer to the following tables for specific information regarding compliance to the standard.

DSK-1031, as a medical electrical equipment, needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

- Portable and mobile RF communications equipments can affect the device.
- The use of accessories other than those specified in this manual may result in increased emissions or decreased immunity of the device.
- DSK-1031 should not be used adjacent to or stacked with other equipment.

Table 201. Guidance and manufacturer's declaration - electromagnetic emissions

DSK-1031 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1031 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR11	Group 1	DSK-1031 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 202. Guidance and manufacturer's declaration - electromagnetic immunity

DSK-1031 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1031 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of DSK-1031 requires continued operation during power mains interruptions, it is recommended that DSK-1031 is powered from an uninterruptible power supply or a battery.
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Table 204. Guidance and manufacturer's declaration - electromagnetic immunity

DSK-1031 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1031 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of DSK-1031, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = [3.5/V1]\sqrt{P}$ <p>$d = [3.5/E1]\sqrt{P}$, 80MHz-800MHz $d = [7/E1]\sqrt{P}$, 800MHz-2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,¹ should be less than the compliance level in each frequency range.² Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m	
NOTE 1:	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
<p>¹ Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which DSK-1031 is used exceeds the applicable RF compliance level above, DSK-1031 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating DSK-1031.</p> <p>² Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 206. Recommended separation distances between portable and mobile RF communications equipment and DSK-1031

DSK-1031 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of DSK-1031 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and DSK-1031 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz to 80 MHz, $d=[3.5/V1]\sqrt{P}$	80 MHz to 800 MHz, $d=[3.5/E1]\sqrt{P}$	800 MHz to 2.5 GHz, $d=[7/V1]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1:	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

CE 0123

NIHON SEIMITSU SOKKI CO., LTD.

2508-13 Nakago Shibukawa Gunma 377-0293 Japan

© Copyright Nihon Seimitsu Sokki Co., Ltd., 2014.

® Registered Trade Mark.

Nissei Healthcare (UK) Ltd. New Barn Lane Henfield West Sussex BN5 9SJ UK

E473/1408/01

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Produsul companiei Nihon Seimitsu Sokki Co., Ltd.